

User Solution: Alphachron™ He Extraction/Measurement Instrument


Copyright © Neo Vista System Integrators Pty Ltd 2004 – All Rights Reserved.

Document Authors: Laurent Raynal (NVSI), Christopher G. Relf (NVSI) and Brent McInnes (CSIRO) Document Release Date: 9th July, 2004

Uranium-helium geochronology is a highly sensitive and cost-effective method of radiometric age dating that can be used to determine the thermal history of the Earth's crust. Helium (^4He), in addition to the radiogenic isotopes of lead (atomic symbol: Pb), is a natural fission product of the uranium (U) and thorium (Th) decay series. Similar to the more commonly used U-Pb geochronology system, age relationships can be determined by measuring the concentrations of both the parent (U and Th) and daughter (^4He) elements in minerals. Australia's CSIRO (Commonwealth Scientific & Industrial Research Organization) is developing applications based on the U-Pb and (U+Th)-He decay schemes to quantitatively determine the thermal histories of mineral belts and petroleum basins, data that is fundamental in the exploration for deposits of minerals, oil and gas.

The Solution

Neo Vista System Integrators Pty Ltd (NVSI) was approached by CSIRO to create a software solution. The task of the software was to tie several hardware technologies together in order to help develop R&D laboratory experimentation into a commercially robust prototype. The result is a professional, industry-focused solution called the *Alphachron™ He Extraction/Measurement Instrument*. Alphachron™ is an automated, integrated and compact turnkey system, designed for the extraction and measurement of gases from mineral samples. CSIRO and its collaborator, Patterson Instruments Pty Ltd, deliver and commission Alphachron™ instruments worldwide. The technologies tied together by the NVSI solution include a Balzers Prisma QMS200 vacuum pump (RS-232), Pfeiffer vacuum heads (TCP/IP), a Eurotherm temperature controller (RS-232), a Coherent laser system controller (GPIB), a Newport motion controller (RS-232), a



SONY camera (through a National Instruments PCI-1408) and numerous vacuum valves (National Instruments PCI- 6503). Two PCs are used to control the system – one in the light-tight room directly controlling the system, and a second in the technician's room for process definition and control, and system monitoring.

When controlling and monitoring sample irradiation, the software solution can automatically detect and correct for inconsistent sample orientations, by determining two reference points, and calculating the relative coordinates of the 25 specimens. Once irradiated with a compact, novel and inexpensive CSIRO-developed diode laser module, the samples begin to glow as they outgas.

